

Non-Absorbable Surgical Suture

EGYTAPE

- Polyester (PET) Tape
- Synthetic
- Multifilament
- Undyed



STERILE | EO

1) Description

- Sterile synthetic non-absorbable polyester tape which is composed of poly ethylene terephthalate (PET)
- Braided
- Available undyed (white)
- Complies with all requirements of USP and EP
- Sterilized using Ethylene Oxide (EO)

2) Indications

Special suture Indicated for the treatment of cervical incompetence during pregnancy–ligation of cervix according to the McDonald technique (purse string ligature). It can also be used in soft tissue approximation or ligation like in orthopedic procedure, as retracting and hanging device as a sling.

Surgery	Needle size
GYNAECOLOGICAL	From 48mm to 65mm

3) Performance

Elicits a minimal initial inflammatory tissue reaction which is followed by gradual encapsulation of the suture by fibrous connective tissue. Not absorbed, nor degraded. Implantation studies in animals show no meaningful decline in polyester suture strength over time.

4) State of Product & Way of Supply

- Sterile
- Strands in size 5 mm
- Half dozen box



Needle			Thread			Code
Needle Description	Shape	length/ Symbol	Color	USP	length	
Without Needle			White	6(5mm width)	1x30 cm	VT8600B
			White	6(5mm width)	1x75cm	MT8600B
Double TP 1/2C	48 mm		White	6(5mm width)	1x30 cm	VT8674B
Double BL 1/2C	65 mm		White	6(5mm width)	1x30 cm	VT8675B BL

